RAD-X

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Cihaz Adı / Device Name: RAD-X Radiological Picture Archiving and Communication System (PACS)

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RAD-X is a flexible archiving, image analysis and reporting system developed to meet the radiology workflow needs of hospitals and hospital groups of all sizes. RADX which offers advanced tools such as 3D reconstruction in customized integrated modules for different modalities such as magnetic resonance (MR), computed tomography (CT), x-ray (DX), mammography (MG), angiography (XA), positron emission tomography (PET), ultrasound (US), MX, RTPLAN, NM, OT, ECG, SR, CR, KO, OPT, PX, PT helps radiologists and clinicians to make the most accurate diagnosis in the fastest time. RAD-X server components save staff time by automating all radiology workflows such as repetitive automatic transmissions within the hospital, reporting distribution, report prioritization, worklist distribution.

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Information/Note: Contains information, advice and tips for better understanding of the instructions given during system operation.

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
Caution: Indicates a hazardous situation that requires professional medical attention, which may result in death, possible injury, or serious illness of the patient.

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RAD-X is a modular system to meet all the needs of radiology departments in hospitals, from image archiving to reporting, from workflow management to directional action research. The RAD-X platform, whose basic philosophy is high integrability, can be used as a complementary system with existing solutions in hospitals as well as providing a total solution.

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DICOM compatible imaging devices are diagnostic imaging devices available in hospitals such as Magnetic Resonance Imaging (MRI), Computed Tomography (CT), Ultrasonography (US), Digital X-Ray (DX), Angiography (XA), Digital Mammography (MX). These devices transfer the images they have obtained to the RAD-X Archive Server or any PACS server in accordance with DICOM (Digital Imaging and Communication in Medicine) standards. DICOM (Digital Imaging and Communications in Medicine) is a digital data format standard developed for the storage, display and analysis of 2D and 3D scientific data received from medical imaging devices.

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Standard development studies are carried out by the "Medical Imaging & Technology Alliance" department within The Association of Electrical and Medical Imaging Equipments Manufacturers headquartered in the United States. Dicom official website is at <https://www.dicomstandard.org/>

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The Dicom standard consists of 17 sub-topics in total, listed below.

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- Part 1: Introduction and Overview
- Part 2: Conformance
- Part 3: Information Object Definitions
- Part 4: Service Class Specifications
- Part 5: Data Structures and Encoding
- Part 6: Data Dictionary
- Part 7: Message Exchange
- Part 8: Network Communication Support for Message Exchange
- Part 10: Media Storage and File Format for Media Interchange
- Part 11: Media Storage Application Profiles
- Part 12: Media Formats and Physical Media for Media Interchange
- Part 14: Grayscale Standard Display Function
- Part 15: Security and System Management Profiles

- Part 16: Content Mapping Resource
- Part 17: Explanatory Information
- Part 18: Web Services
- Part 19: Application Hosting
- Part 20: Imaging reports using HL7 Clinical Document Architecture
- Part 21: Transformations between DICOM and other Representations
- Part 22: Real-Time Communication

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Parts of the DICOM Standard

- Part 1: Introduction and Overview
- Part 2: Conformance
- Part 3: Information Object Definitions
- Part 4: Service Class Specifications
- Part 5: Data Structures and Encoding
- Part 6: Data Dictionary
- Part 7: Message Exchange
- Part 8: Network Communication Support for Message Exchange
- Part 10: Media Storage and File Format for Media Interchange
- Part 11: Media Storage Application Profiles
- Part 12: Media Formats and Physical Media for Media Interchange
- Part 14: Grayscale Standard Display Function
- Part 15: Security and System Management Profiles
- Part 16: Content Mapping Resource
- Part 17: Explanatory Information
- Part 18: Web Services
- Part 19: Application Hosting
- Part 20: Imaging reports using HL7 Clinical Document Architecture
- Part 21: Transformations between DICOM and other Representations
- Part 22: Real-Time Communication

DICOM compatible devices transfer the image files created in accordance with the 10th part of the DICOM standard to the archive server in accordance with the network communication and message sending rules described in the 8th part of the standard.

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Endoscopy, colonoscopy, document scanner, etc. devices are devices that obtain images that are not in DICOM standard. The images obtained by these devices can be in many formats such as Bitmap, JPEG, AVI, PDF, etc. Since these devices are not DICOM compatible, there is no patient and examination data in the images. Therefore, these images must be matched with certain patients and examinations within the framework of certain rules on the RAD-X archive server. The RAD-X archive server provides manual matching as well as automatic rules definition for this matching.

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They are systems that query the examinations on DICOM archive servers and take the images of the desired examination in DICOM format. In this group, workstations that offer advanced image processing tools come first. For example, a system for spectroscopy analysis on MRI images provides the opportunity for radiologists to work by taking the images to be examined from archive systems. Since these systems generally do not have long-term archiving capacity, the analyzed images are deleted after a certain period of time and downloaded again when necessary. Communication standards in DICOM Part 8 are used for communication with archive servers.

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RAD-X archive server reads and stores metadata of patients, examinations and images in DICOM images from imaging devices and other sources. The files are also put on the file system and the address reference information is kept in the database. RAD-X archive server provides the following service features mentioned in DICOM standards.

- C-Store SCP: Service that receives images sent from devices
- C-Find SCP: Service that enables devices and workstations to query the examinations in the system
- C-Move SCP: Service that fulfills requests of devices and workstations to send certain images to another source
- C-Echo SCP: Service that allows other DICOM compatible systems to query whether the RAD-X server is running
- DICOM MWL SCP: System that allows imaging devices to receive patient and examination information electronically
- DICOM Print SCU: System that allows images to be sent to DICOM printers that support the DICOM Print protocol

RAD-X archive server can use MSSQL, MySQL, PostgreSQL and Oracle systems as database.

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HIMS are systems where the radiology workflow starts and ends. When a radiology examination is requested for any patient, this process is done through HIMS. HIMS transmits radiology requests by integrating with the RAD-X reporting system and provides feedback to the relevant physician by taking back the reports prepared.

RAD-X Raporlama Sunucusu / RAD-X Reporting Server

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RAD-X Reporting server acts as a bridge between HIMS and PACS systems. Its main function is to receive the request information created in HIMS and then to establish a relationship between these requests and the images in PACS systems. Thus, any radiologist will be able to write a report on the right request by looking at the right image. The RAD-X reporting server contains tools that allow radiologists to generate reports in the best and fastest way. For example, by offering different report templates for different examinations, radiologists can write reports quickly. Thanks to the structured reporting methods, which have become increasingly important recently, radiologists can prepare reports faster and without skipping any point by choosing the appropriate options for the questions directed to them instead of creating text reports for any examination. RAD-X reporting server provides this infrastructure to users thanks to its expandable structured report options. In addition to the text-based reporting option, the voice reporting option is also available in RAD-X. Thanks to this application, radiologists can create their reports simply by speaking. These voice recordings are then transcribed by the secretaries and made suitable for sending to HIMS.

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RAD-X Reporting server supports HL7 (Health Level 7) message standard for communication with HIMS. HL7 refers to Health Level Seven Inc., as well as some standards developed in this context (eg HL7 v2.x, v3.0, HL7 RIM, etc.), which are entirely voluntary and intended to establish international health standards. It is developed by a team headquartered in Ann Arbor, Michigan, USA, by the Standards Developing Organization (SDO) authorized by the American National Standards Institute (ANSI). Established in 1987 to develop a standard for Hospital Information Management Systems (HIMS), HL7 is the only standard preferred by many institutions and organizations in the transmission of clinical data today. HL7 is dedicated to developing a comprehensive framework for the exchange, association, sharing and use of electronic health data. Using this standard, the RAD-X Archive server receives radiological request information from HIMS (HL7 ORM) and transfers the generated reports to HIMS (HL7 ORU).

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RAD-X Viewer Platform is a web-based system that allows querying of examinations on all kinds of PACS servers, especially RAD-X Archive Server, and accessing related images. The RAD-X Viewer Platform differs from other systems in that it provides access to different image sources. Querying the audits is performed on the server side with the C-Find services described in DICOM Chapter 8. Query parameters entered by users from the web interface are converted into DICOM format and transmitted to PACS servers that offer C-Find service. Query responses from these systems are displayed in the appropriate format on the web-based viewer. Image accesses are provided by the Web Access to DICOM Persistent Objects (WADO) services described in DICOM Chapter 18. Thus, images from image sources are transferred directly to the client screen. RAD-X Viewer Platform supports all standard image formats such as MR, CT, DX, MX, US, XA, RTPLAN, NM, OT, ECG, MG, SR, CR, KO, OPT, PX, PT. It offers tools that allow all standard measurement functions such as length, angle, density, area, as well as all standard image processing functions such as windowing, 3D positioning, magnification, geometric transformations on these images.

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An important difference of the RAD-X Viewer Platform from other systems is that it contains different modules (RAD-X Viewer plug-ins) for different image types. The most suitable plug-in for the image type in the archive is automatically loaded and presented to the user.

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The RAD-X Viewer Platform is a class IIb medical device within the scope of the Medical Device Regulation and can only be used **by authorized users**.

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RAD-X can be accessed via a web browser using the secure access address determined by the institution.

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The following requirements exist for desktop computers;

- Display with at least 1366x768 resolution (Recommended 1920x1080 and above)
- At least 4GB of RAM (8GB or more recommended for multislice images such as CT and MR)
- At least 2 GHz processor
- At least 10GB free disk space
- Chrome, Firefox or Edge browsers
- Stereo microphone (for sound recording)

The following requirements exist for mobile devices;

- iOS 12 or higher operating system
- Android 9.0 or higher operating system

Access must be provided over a secure connection for voice recording and 3D imaging. For voice recording, access to the microphone must be granted through the browser.

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RAM : Minimum 8 GB

Disk : Minimum 1.5 TB

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[Minimum system requirements for the server running the application services:](#)

[Operating System: Windows Server 2016 or later](#)

[CPU: Minimum 8 Cores](#)

[RAM: Minimum 8 GB](#)

[Disk: Minimum 1.5 TB](#)

[Network Connection: 1 Gigabit Ethernet](#)

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[The following services run on the application server:](#)

- [Web applications](#)
- [API services](#)
- [Authentication services](#)
- [System integration services](#)

[SSD or NVMe disk infrastructure is recommended to improve system performance.](#)

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RAM : Minimum 8 GB

Disk : Minimum 550 GB

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[Minimum system requirements for the server running the database services:](#)

[Operating System: Ubuntu Server 22.04 LTS or later](#)

[CPU: Minimum 6 Cores](#)

[RAM: Minimum 8 GB](#)

[Disk: Minimum 550 GB](#)

[Network Connection: 1 Gigabit Ethernet](#)

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The following services run on the database server:

- PostgreSQL or relational database system
- Data processing services
- System logging mechanisms

Using an SSD disk is recommended for high-performance data access.

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The system has architecture that operates over a local area network (LAN).

Minimum network requirements:

Protocol: IPv4

Network Type: Ethernet

Minimum speed: 1 Gbit/s

Recommended speed: 10 Gbit/s

The system can operate over a 1 Gigabit Ethernet connection. However, a 10 Gigabit Ethernet infrastructure is recommended to improve performance in cases of large data transfers and heavy user load. Wireless technologies are not used in the system design.

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Data communication between system components is carried out using the following protocols:

Network Layer: IPv4

Transport Layer: TCP

Application Layer: HTTPS / REST API

Data transfers are carried out via TCP-based communication.

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The following security mechanisms are used to ensure the security of data transmitted within the system:

- SSL/TLS encryption
- HTTPS communication protocol
- Server certificate verification
- Authentication services

Data communication between system components is encrypted via the SSL/TLS-based HTTPS protocol. This reduces the following security risks:

- Packet sniffing
- Man-in-the-middle attacks
- Unauthorized data access

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Radiologists review and report radiological images in the light of clinicians' findings and preliminary diagnoses. Radiologists are divided into different branches, especially in large-scale hospitals. For example, while different radiologists perform MRI and CT scans in some hospitals, larger hospitals have radiologists specializing in different anatomical regions such as brain, thorax, abdomen, and joints. The tools used by radiologists in each examination may differ. That's why the RAD-X Imaging Platform offers a customizable array of screens and tools for different inspections. Radiologists also prepare reports for clinicians' use using the reporting tools provided by the RAD-X reporting system.

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It is the term used for doctors who frequently apply to image-based diagnostic methods such as neurology, orthopedics, chest diseases, emergency, cardiology. Interaction of clinicians with the radiology system begins at the stage of requesting examination. After the images of the patients are taken and radiological reports are created if necessary, clinicians try to make a definitive diagnosis by accessing these data and evaluating them together with other findings. Clinicians can access images of their patients, reports and if any radiologists' markings and comments on images containing pathology via the RAD-X Viewer Platform.

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The administrator is the administrative staff authorized to receive reports on the overall system and the technical manager staff who are authorized to take the actions that require managerial level intervention regarding the operation of the system. RAD-X provides independent management screens that only users in this group can access. Operations such as system configuration, user management, component management, process management, and administrative reporting can be performed on these screens.

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
Other users are personnel who are not physicians such as radiology technicians, medical secretaries, data entry personnel, but who use RAD-X systems for purposes such as querying, preparing reports, taking printouts and obtaining information to the extent of their authority. By creating different user profiles in RAD-X, it is possible to authorize these profiles separately.

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Username and password required to access RAD-X are given to users by system administrators. Passwords must be a minimum of 8 characters; must consist of at least one lowercase character, one uppercase character, one number, and one non-numeric or non-letter character. Password change is mandatory at intervals determined according to your institution's data policies. The RAD-X access screen is shown below.

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

The RAD-X web interface consists of the components shown in the figure below. See Detailed Interface Components 3234 for details of components .

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

B I U S " U L </> H1 H2 Normal 1.5em

K Templates Save Template

During the MRI examination of your patient's Upper Abdomen,

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
For more detailed information, see Annex-2: Display Toolbar.

⑥ Menü Çubuğu / Menu Bar

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
Symbol	Command	Explanation
	3D Distance Measurement	It is used to measure the distance between two points (in MPR mode with other projections) on a 3D image.
RG	Region Growing	It is used to delete within a spherical region determined on the image by region growing method.
FHC	Free Hand Cut	It is used for volumetric deletion of intersections along the image axis of an area determined by freehand drawing.

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For more detailed information, see Annex-3: Menu Bar.

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
The system will place the components on the screen in the most appropriate way according to the user type and screen resolution. If the user wishes, user can choose different screen layout options from the menu that opens with the **menu** (≡) button at the top left and is shown in the figure below.

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To exit RAD-X safely, you can exit the system by pressing the **“Exit”** button in the menu bar.

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
Frequently used search criteria are presented to the user in the study search window as shown in the figure below. After the desired criteria are entered, the search results are listed by pressing the **magnifying glass** (🔍) button. The **Clear button** is used to clear the search after adding criteria.

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

The criteria can be edited at the request of the user. After pressing the **arrow** (▼) button under the magnifying glass, the desired criteria can be added or removed as seen in the figure below.

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Maximum 4 criteria can be added at the same time. Use ⊕ to **add** , ⊖ to **subtract** .

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As seen in Figure 12, detailed searches can be made with the Additional Criteria List shown in Figure 13, which appears with the arrow sign under the magnifying glass for patient searches with more than 4 criteria. After entering the information, the SEARCH button should be pressed.

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In order to repeat searches with certain criteria, you can name your search by pressing the **"Save last search"** button, then you can find it by clicking the vertical ellipsis (**⋮**) in the search section.

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After the desired criteria are added or removed, the system always saves the search criteria as these criteria for the user in the system when the **"Save search fields"** button is pressed.

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
PACS Driven mode is designed for mismatched and urgent examinations. Reporting cannot be performed on examinations opened in PACS Driven mode.

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RAD-X offers training folder options for easy labeling of exams and easy access when required. By clicking the **file** (■) button in the search window, the training folder menu is displayed as in the figure below. From this menu, operations such as creating/deleting new folders, sharing the created folders, assigning open examinations to the desired folders are performed. To access the recordings in any training folder, it is sufficient to click on the folder.

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
RAD-X offers saved worklist options to allow users to work quickly. By pressing the **vertical ellipsis** (⋮) button in the search window, the worklist options are displayed as in the figure below. Clicking on the desired worklist option will perform the search.

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Worklist options are predefined according to user type. Users can expand these job listings if desired. For detailed information Check out the Worklist.

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

Wrong e-mail or password: The user must use the username and password defined in the system.

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
Error Code 500: This occurs when the server encounters an unexpected condition that prevents it from fulfilling the request. This issue cannot be resolved by the user and requires server-side intervention by a software developer or system administrator.

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
Error Code 501: It occurs that the server does not support or has not implemented the functionality required to fulfill the request. In other words, the server has received the request but does not know how to process it, or the relevant feature has not yet been activated. This error cannot be resolved by the user and requires server-side intervention by a software developer or system administrator.

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Error Code 502: This indicates that an invalid or faulty response is received from the server due to the expiration of the software license. This error cannot be resolved by the user and requires server-side intervention by a software developer or system administrator.

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Code 400 error: Indicates that the server cannot understand or process the request. This is usually caused by an incorrect, incomplete, or improperly formatted request sent by the client. The user cannot resolve this error. Intervention on the server side by a software developer or system administrator is required.

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In case you need maintenance and assistance regarding the application, you can contact us by using one of the communication tools given in the contact information section. Support requests will be answered within 24 hours, and priority management will be completed and decided within 48 hours. Support will be given by Simplex Information Technologies Inc. support team. The software team will take charge in unresolved issues.

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

For detailed information about the buttons used in the Examination and Work List, please refer to Annex-1: Buttons.

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If users are authorized to open images, patient images are opened in the viewing window shown in the figure below by double-clicking a record in the examination list or by pressing the Open Image button ().

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
Examinations in all hospitals belonging to the patient are listed in the patient history window. The most suitable layout for the analyzed study is automatically selected. Clicking the "HP" button on the toolbar enables quick switching between predefined trigger-specific protocols. In order to display any series in a window, one of the following methods can be followed with the stamp images seen in the Figure 32 Image Information Window.

- Stamp image belonging to the series can be dragged and dropped to the desired window.
- After clicking on the desired window, the stamp image to be opened can be clicked once.
- The related cell can be clicked with the left mouse button from the graphical menu that appears when the serial stamp image is right clicked. When a cell in this menu is clicked with the right mouse button, all series are automatically distributed to the windows so that the selected series is to the first window by automatically changing the sequence.

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When any image window is selected, the most suitable module is automatically selected according to the image type in that window

and the buttons of that module are presented to the user in the toolbar.

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RAD-X contains the following modules.

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
RAD-X displays each series in an examination in a viewport. It is possible to open up to 64 windows on the screen. Image-specific labels are displayed in the four corners of each window. These tags can be set modality specific. The **ANN** button in the toolbar is used to open, close or partially close the tags.

The scale is shown next to the display window. If the scale value of a pixel is recorded in the picture, this scale is shown in mm as an example. Otherwise, it is displayed in pixels. To turn the scale on/off, the scale ($\leftrightarrow$) button on the toolbar can be pressed.

With the buttons at the bottom of the display window, it is possible to play the series in cine mode, switch to the previous/next series, download the displayed image or copy it to the clipboard. In addition, you can see

the original, axial, sagittal, coronal sections of the image, put it in MIP mode, perform a 4-section study and examine the report.

There are three different settings (MinIP, AvgIP, MIP) in MIP mode and these settings are indicated at the bottom right of the display window. In addition, the numbers for the calculation of thickness (thickness of the image x the number of images) are located at the bottom right. A maximum of 8 slices can be used in MIP mode. The **thickness** ($\frac{1}{n}$) button is used to change the number of sections 1-8.

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- › The RAD-X Imager Platform is a class IIb medical device within the scope of the Medical Device Regulation and may only be used by authorised users.
- › PACS Driven modu eşleşmeyen ve acil olarak incelenmesi gereken tetkikler için tasarlanmıştır. PACS Driven modu açılan tetkikler üzerinde raporlama işlemi yapılamamaktadır.
- › PACS Driven mode is designed for tests that do not match and need to be examined urgently. Reporting cannot be made on examinations for which PACS Driven mode is switched on.
- › Radyologların ses kaydı yaparken kısaltma kullanmaları, diktasyonda hatalara yol açabilir.
- › The use of abbreviations by radiologists during voice recording may lead to errors in dictation.
- › Ses kaydının başarı ile tamamlanması için tarayıcı üzerinde izin verilmesi gerekmektedir. Eğer aktarım bitmeden tarayıcı kapatılırsa, ses kaydı gönderilemez.

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For detailed information on the buttons (toolbar) used in the Display Window, please see Annex-1: Buttons and Annex-2: Display Toolbar.

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
On the display screen, when one of the patient's old examinations is opened, "prior" is written in the upper left corner, "current" is written when a new examination is opened, and “post” is written if a new examination is opened from an old examination.

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These functions, selected according to the type of examination, can be assigned to different buttons of the mouse and are functions that make changes on the image by moving the mouse while holding the mouse button on the picture. The mouse buttons that can be used and their codes are shown in the picture below. The desired function is selected from the toolbar by clicking the mouse button to be assigned. Which function is assigned to which button is indicated by the number codes on the buttons.

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
In RAD-X, the VRTViewer option appears at the bottom right of the series that are available for 3D viewing.

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
After the serial is converted to volume, the viewer switches to 3D review mode.

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In 3D display mode, position can be adjusted on the right side of the window and thickness can be adjusted on the left side. Coloring options are accessed from the TF window in the toolbar. By selecting VTR mode from the HP window, MPR mode is entered. All windows can be synchronized with 3D Crosshair in MPR mode.

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
For detailed information on the buttons (toolbar) used in the 3D Display Window, see Annex-2: Display Toolbar.

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RAD-X allows the creation of one or more examination reports based on the opened image. The reporting window consists of 4 parts as shown in the figure below.

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- › A serious incident related to our device includes:
 - Death of a patient or user
 - Serious deterioration in the state of health of a patient (e.g. permanent impairment, life-threatening condition, hospitalization or prolongation of hospitalization)
 - A serious public health threat
- › In the event of such an incident, it must be reported immediately to the manufacturer, the Competent Authority of the European Union Member State in which you reside, and the notified body. The manufacturer information and contact details are provided below:
 - **Simplex Bilgi Teknolojileri A.Ş. /**
BAHÇELİEVLER MAH. DOKTOR CEVDET
KARA CAD. ANKARA ÜNİVERSİTESİ
TEKNOKENT NO: 35 C İÇ KAPI NO: 112
GÖLBAŞI/ ANKARA
 - Phone: +90 312 485 65 05
 - Email: info@simplexbt.com

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


For detailed information on the buttons used in the Reporting Window, see Annex-1: Buttons.

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From this section, the examinations to be reported are selected. All connected examinations of the opened image are listed in this section with their SUT codes and status information. If reports are to be written to more than one examination at the same time, these examinations are selected from the selection box on the left.

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Voice recording can be activated by pressing the  button. Press the same button again to complete the voice recording. Each recording track is listed in this section. Desired recording tracks can be played back, deleted or overwritten. Once the radiologists approve the recording pieces they have created, all the pieces are combined and transmitted to the reporter in the form of a single sound recording.

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The use of abbreviations by radiologists when recording voice can lead to errors in dictation.

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Permission must be granted on the browser for the voice recording to be completed successfully. If the browser is closed before the transfer is finished, the voice recording cannot be sent.

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With the editor in this section, text reports in HTML format can be prepared. Saved templates can be loaded from the Templates menu or a written text can be saved as a template. Report templates include findings, method, comparison, and results headings. A report is created by entering the relevant text under each heading.

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Because reports are sent to the Hospital Information System after approval, radiologists must check before confirming.

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
The Examination Status Codes specified in the table below are described in diagram form in the **Reporting Scheme - Buttons and Actors** and **Reporting Workflow Sections**.

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


Figure 40 Reporting Schemas for Assistances

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The Assistant Radiologist can perform all the operations performed by the radiologist, except for approving the report (Code: 180).

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Query	Explanation	Code(s)	Criterion	Persons
To be reported	Examinations that are ready to be reported but not assigned a radiologist.	100		Radiologist
To be reported by me	They are ready to be reported and assigned to a radiologist.	110	Assigned Radiologist	Radiologist
Reported	All approved reports.	180		Radiologist, Rapporteur
Reported by me	Approved reports specific to each radiologist.	180	Assigned Radiologist	Radiologist
To be reported urgently	These are examinations that need to be reported urgently and that have not been assigned a radiologist.	100	Priority code should be checked (20-Emergency)	Radiologist
To be confirmed	A list of all reports to be approved. Approval of assigned radiologists is awaited. (to be confirmed by me)	160		Radiologist
To be confirmed by me	These are the reports that the assigned radiologist must confirm.	160, 170, 175	Assigned Radiologist	Radiologist
Voice recorded by me	These are the voice recordings that the assigned radiologist sent to the system before the report was approved.	130, 140	Assigned Radiologist	Radiologist
Additional procedure to be performed	Additional procedure is the desired patient list.	162, 163		Technician
Additional procedure to be confirmed	It is the patient list that the technician will approve the procedure after the additional procedure is done.	164, 165		Technician
New voice record	It is the report list that needs to be written, the radiologist made a voice recording.	130,140		Rapporteur
To be dictated by me	These are the reports that the assigned rapporteur should write.	140	Assigned Rapporteur	Rapporteur
Dictated by me	They are written reports specific to each rapporteur.	160	Assigned Rapporteur	Rapporteur

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- Serious deterioration in the state of health of a patient (e.g. permanent impairment, life-threatening condition, hospitalization or prolongation of hospitalization)
 - A serious public health threat
- In the event of such an incident, it must be reported immediately to the manufacturer, the Competent Authority of the European Union Member State in which you reside, and the notified body. The manufacturer information and contact details are provided below:
- **Simplex Bilgi Teknolojileri A.Ş /** BAHÇELİEVLER MAH. DOKTOR CEVDET KARA CAD. ANKARA
ÜNİVERSİTESİ TEKNOKENT NO: 35 C İÇ KAPI NO: 112 GÖLBAŞI/ ANKARA
 - Phone: +90 312 485 65 05
 - Email: info@simplexbt.com

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Symbol	Command	Explanation
	Details	Displays the patient's information.
	Send DICOM	It allows the patient's images to be sent to a different server or CD Robot.
	Upload Image	It is used to upload images of the patient from other hospitals.
	Open Image	Displays images of the patient.
	Specify Columns to Display	It is used to change the columns shown in the survey list as desired.
	Download CD Image (ISO)	It is used to download examinations in ISO format.
	Refuse Recording to Text	It is used to give up the reports that rapporteurs have received.
	Open Approved Report	It is used to open the report approved by the radiologist.
	Change Priority	It is used to change the priority of the patient. (eg VIP, Emergency, Normal)
	Cancel Reporting	It is used to give up the reports that radiologists have received.
	Open Report	It is used to see the reports that have been written or are being written.
	Download Table in Excel Format	It is used to download the listed patient examination list as a table in Excel format.

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Symbol	Command	Explanation
	Preload	Allows pre-downloading images of a test before it is opened. In this way, the downloaded images are loaded faster when the examination is opened. It is recommended to use this command for a maximum of 3 exams, as preloading too many exams can cause performance problems. After opening a previously uploaded examination, the downloaded images are deleted when a new examination is opened.
	Assign to Radiologist	It is used by the specialist radiologist to assign to a radiologist without opening an examination. Authorized users can re-pool an examination from a 110 state assigned to a radiologist with this command and draw it to the 100 state.

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Symbol	Command	Explanation
+	Add file	It is used to upload a file into the examination. Each uploaded file appears in stamp images as a new series. Files in PDF, PNG, MP4, JPG and BMP formats can be viewed directly in the viewport. For other file types, the user is given the option to download to the computer.
i	Patient Information	Shows information about the patient and the selected request.
💬	Patient Messages	It allows users to enter additional information about the patient. The entered information is only displayed by the user who entered the message or by the users addressed on the message screen.
≡	Reports	Shows reports of requests associated with the selected study. Reports can be copied from the pop-up window.
🎨	Measurements	Shows the list of measurements made on the images of the examination. Relevant measurements can be deleted from this list or a quick transition to the image containing the measurement can be made.
🔗	Save Measurements Generally	It is used to record the measurements made in the examination images in a way that everyone can see. Measurements made in this way are automatically loaded when the survey is opened by any user.

	Save Measurements Specifically	It is used to record the measurements made in the examination images in a way that only the active user can see.
	Print Image with Measurement	It is used for whole printing of the images measured in the examination.
	Reload Measurements	Reloads saved measurements.
	Delete All Measurements	It is used to delete all measurements made on images in the examination.

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Symbol	Command	Explanation
	Print with Subtitle	It is used to print the approved report by adding special subtitles under it.
	Send Back to Assistant	It is used for the rejection of the report pre-approved by the assistant by the radiologist and for the assistant to rewrite the report.
	Detail Information Screen	For statistical purposes, it is used to select the report or finding category defined by the institution.
	E-Signature	It is used for the radiologist to electronically sign the report. (For e-signature, electronic signature device must be connected to USB port and electronic signature application must be installed.)
	Patient Admission for Additional Procedure	It is used by the technician to indicate that the patient is taken to additional procedure.
	Additional Procedure Cancel	It is used for the radiologist to cancel the additional procedure request.
	Request Additional Procedure	It is used by the radiologist to request additional procedure.
	Additional Procedure Confirmation	It is used by the technician to indicate that the additional procedure has been completed after the additional procedure has been performed and the images have been sent to PACS.
	Submit Additional Report	It is used to send the written report in addition to the approved report.
	Submit Additional Voice Recording	It is used to send the additional voice recording made by the radiologist.
	Record Additional Voice	It is used by the radiologist when additional voice recording is desired.
	Cancel Additional Voice Recording	Used to exit the additional voice recording mode.

	Image Feedback	It is used by the radiologist who wrote the report to provide feedback on image quality.
	Image Quality-Control	It is used for evaluation of the image by an independent quality-control specialist assigned for quality-control.
	Confirm and Send to HIMS	It is used to approve the report written by the assistant radiologist/rapporteur by the radiologist.
	Confirm and Send to HIMS	It is used by the radiologist to confirm the report.
	Edit Approved Report*	It is used to edit the report approved by the radiologist.
	Select Expert to Approve	It is used to select the radiologist who will approve the report written by the assistant radiologist.
	Previous Reports	It is used to see previous reports of the patient.
	Send to Radiologist	It is used to send the report written by the rapporteur to the radiologist.
	Report Quality-Control	It is used for evaluation of the report by an independent quality-control specialist assigned for quality-control.
	Cancel Reporting	It is used by the radiologist to remove the assigned exam from the work list without reporting.
	Start Reporting	It is used to add the examination to the radiologist's worklist.
	Send Back to Rapporteur	It is used to reject the report written by the reporter by the radiologist and to rewrite the report by the reporter.
	Undo Voice Recording	It is used to withdraw the voice recording sent by the radiologist from the work pool of the reporters.
	Send a Voice Record	It is used to send the voice recording recorded by the radiologist to the work pool of the reporters to be converted into text.
	Save Draft	It is used to save the report to continue later without moving on to the next job step.
	Undo All Operations*	It is used to undo all the operations performed on the patient.
	Cancel All Report Operations*	It is used to cancel the report approved by the radiologist. Voice recordings and written reports are deleted with this button.
	Submit for Expert Approval	It is used to send the report written by the assistant radiologist to the radiologist's approval.
	Print	It is used to output the report of the images.
	Stop Writing	It is used by the reporter to remove the assigned examination from the work list without reporting.
	Start Writing	It is used to add the examination to the secretary's work list.

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Symbol	Command	Explanation
FPS	Number of segments played per second	It is used to change the playback speed of multi-section series in CINE mode.
ANN.	Texts on the image	It is used to adjust the detail level of texts containing patient and examination information on the image (eg All, Partial, Off).
TB.	Toolbar	It is used to set the display level of the functions used for the image (eg All, Partial, Off).
HP	Hanging Protocol	It is used to implement hanging protocols defined by examination type.
WL	Window Level	It is used to apply window settings (eg Abdomen, Bone, Brain) created specifically for examination to the selected CT image.
SYNC	Mammography sync mode	It is used to set the synchronization status of the brightness and magnification of mammography images.
CB	Contrast and Brightness	It is used to apply the contrast and brightness settings created specifically for the examinations to the selected MR image.
TF		It is used to select the desired tissue modes for 3D imaging. Coloring
MODE	Display Mode	It is used to select modes specific to 3D display.

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Symbol	Command	Explanation
	Contrast/Brightness	It is used to make Right-Left Contrast and Up-Down Brightness adjustments. It is used to apply brightness and contrast to the entire sequence in case of images with different window settings in the sequence, such as MR.
	Window Level	It is used to make Right-Left window width and Up-Down window level settings. Changes the window setting of the entire sequence in the case of images with the same window settings in the sequence, such as CT. It only changes the window of an image in case of images with different window settings in the sequence, such as MR.
	Scroll	It is used to switch between slices in multislice images (eg CT, MR).
	Zoom	It is used to enlarge or shrink the image size.
	Pan	It is used to move the image.

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Symbol	Command	Explanation
	3D Crosshair	If the selected point in the active image is among other sequences open on the screen, these sequences are used to mark the position in that segment by wrapping them in the segment containing the point. If the selected point is not contained within $\pm 2\text{mm}$ of other sequences, the sequence is wrapped in the section closest to this point.
	Thickness	It is used to manually adjust the thickness in MPR mode.

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Symbol	Command	Explanation
	Free Rotate	It is used to rotate the image.

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Symbol	Command	Explanation
	3D Crosshair	If the selected point in the active image is among other sequences open on the screen, these sequences are used to mark the position in that segment by wrapping them in the segment containing the point. If the selected point is not contained within $\pm 2\text{mm}$ of other sequences, the sequence is wrapped in the section closest to this point.
	Thickness	It is used to manually adjust the thickness in MPR mode.
	3D Rotation	It is used to rotate the image in 3D.
	Add CPR Point	It is used to specify points on the curved axis for the Curved MPR mode.

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Symbol	Command	Explanation
	Clear Last Measurement	It is used to delete the last measurement made on the image.
	Clear All Measurement	It is used to delete all measurements made on the image.
	Synchronization	It is used in 3D examinations such as MRI, CT to turn on and off the synchronization of series in the same plan.
	3D	It is used to establish a connection to the native 3D server application.
	Automatic HP on/off	It is used to turn on and off the status of the application of automatic hanging protocols suitable for examinations.
	LUT on/off	It is used to turn the LUT representations of images containing LUTs on or off.
	Send Email	It is used to send images via mail.
	Open file	It is used to download DICOM external files that are subsequently uploaded to the examination to the computer.
	Show/hide scale	It is used to turn on and off the ruler on the edge of the image.
	1-1 Zoom	It is used to resize the image to its original size.
	Invert Colors	It is used to invert the colors in the image.
HFLIP	Flip Horizontal	It is used to flip the image horizontally (relative to the y-axis).
VFLIP	Flip Vertical	It is used to flip the image vertically (relative to the x-axis).
	Reference Lines	In 3D examinations such as MRI and CT, it is used to turn on and off the display of reference lines of the active window on other windows.
	Interpolation on/off	The interpolation method used in the enlargement display of images is used to turn on and off.
	Reset Picture	It is used to undo all the operations done and return the image to its original state.
	Print	It is used to output image or images.
	Measuring Handles	It is used to set the display state of the handles so that measurements can be modified more easily.

BT (CT), MR Ek Fonksiyonlar / CT, MR Additional Functions

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Symbol	Command	Explanation
	Split Pictures Window Settings	It allows the window setting made in multi-section series to be set only to the selected image or to the whole series.
	Create Series and Submit to PACS	It is used to create a new DICOM file as shown in a processed and measured image.

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Symbol	Command	Explanation
	Flip 90 Degrees	It is used to rotate the image 90 degrees clockwise.

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Symbol	Command	Explanation
	True Size Zoom	It provides a representation of the image on the screen in its actual physical size. (The PPI value for the display must be entered correctly)

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Symbol	Command	Explanation
	Clear Notes	It is used to delete the note on the image.
MESH	Create SSD Mesh	It is used to create superficial data from volume data with SSD method. (This function must be run on the volume for SSD and SSD-WF displays.)
	Accept Deletion	It is used to delete the area determined by Free Hand Cut.
	Deny Deletion	It is used to cancel the deletion of the area determined by Free Hand Cut.
	Cancel All Deletions	It is used to get back all areas deleted with Free Hand Cut.
S, I	Superior, Inferior	It is used for the top and bottom view of the image respectively.
L, R	Left, Right	It is used for left and right view of the image respectively.
P, A	Posterior, Anterior	It is used for the back and front view of the image respectively.
	Shading On/Off	It is used to turn the shading of the image on or off.
	Cutting Plane On/Off	It is used to activate the 3D cut box.
	Clear CRP Points	It is used to delete added CRP points.
	Recording Starting Point	It is used to select the video starting point.
	Recording End Point	It is used to select the video endpoint.
	Record Video	It is used to record the video after selecting the start and end point.

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Symbol	Command	Explanation
	Lense	It is used to enlarge a certain area of the image.
	Distance Measurement	It is used to measure the distance between two points.
	Cross Distance Measurement	It is used to make two consecutive length measurements.
	Stenosis Measurement	It is used to measure the stenosis ratio over two elliptical fields.
	Add note	It is used to annotate the image.
COBB	COBB Angle Measurement	It is used to measure the angle between two lines.
	Add Volume Reference Point	It is used to mark spatially the same points on a picture of each of the two examinations in order to synchronize 3D examinations such as CT and MR obtained at two different times.
	Point ROI Measurement	It is used to measure ROI pointwise.
	Automatic Contrast	It is used to make automatic contrast according to the specified area.
	Advanced Settings	It is used to open the advanced settings window specific to the examination type.
	Creating Keyboard Shortcuts	It is used to create keyboard shortcuts.

BT (CT), MR Ek Ölçümler / CT, MR Additional Measurements

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Symbol	Command	Explanation
	DDG Measurement	It is used to measure the distance between two points on a line.
	TT-TG Measurement	It is used to measure TT-TG distance.
	ROI Measurement	It is used to measure the mean, min and max HU or intensity value (ROI) in a rectangular area.
	Elliptic ROI Measurement	It is used for ROI measurement in an elliptical field.
	Free ROI Measurement	It is used for ROI measurement in a free-handed area.

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Symbol	Command	Explanation
SHELL	Add Shell Template	It is used to add predefined shell templates to the image at scale.
CUP	Add Cup Template	It is used to add predefined cup templates to the image in a scaled manner.
STEM	Add Stem Template	It is used to add predefined stem templates to the image in a scaled manner.

CAL	Perform Linear Calibration	It is used to calibrate the image by entering the unit of length between two selected points.
^	Angle Measurement	It is used to measure the angle between two lines.
DIF	Length Difference Measurement	It is used to measure the difference in length between two lines.

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Symbol	Command	Explanation
CAL	Perform Linear Calibration	It is used to calibrate the image by entering the unit of length between two selected points.

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Symbol	Command	Explanation
	Patient List	It is used to open the Examination and Worklist.
	Reporting	It is used to open Reporting Window.
	Stamp Pictures	It is used to open the stamp images window of the examination.
	Display Window	It is used to open the Display Window for the examination.
	JPG Download as DICOM or JPEG	It is used to specify whether images are downloaded as DICOM (left) or JPEG (right).
	Examination Download List	It is used to view the list of exams that have already been loaded.
	Önceden Yüklenen Bir Sonraki Tetkiği Aç	It is used to open the preloaded examinations sequentially.
	CD Download List	Used to view the list of exams downloaded in ISO format.

	Examination Submission List	It is used to view the list of patient images sent to a different server or CD Robot.
	Park	When you log in to the application again by saving the information of the open examination, it is used to resume where you left off with the Restore command.
	Restore	It is used to open the parked examination in the same way.
M1	Monitor 1	It is used to open a separate display window for each monitor.
M2	Monitor 2	
M3	Monitor 3	
PPI	PPI	It is used to adjust the physical screen resolution. (Used for the True Size Zoom feature on mammography images).
	Open Help File	Used to open help videos.
SIGNOUT	Exit	It is used to log out of the system.
	Language selection	It is used to switch between Turkish, English and Russian languages.

1	2	3	4
	Cihaz Adı/Modeli Device Name/Model	 1984	
5	6	7	8
	 XXXXXXXX		
9	10	11	
TF.01.LT.01 / Rev.00 --/--/---	Cihaz Versiyon Numarası Device Version No		Simplex Bilgi Teknolojileri A.Ş Adres / Address: BAĞÇELİEVLER MAH. DOKTOR CEVDET KARA CAD. ANKARA ÜNİVERSİTESİ TEKNOKENT NO: 35 C İÇ KAPI NO: 112 GÖLBAŞI/ ANKARA Telefon / Phone: + 90 312 485 65 05 Web / Mail: www.simplexbt.com / info@simplexbt.com

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© Copyright, Simplex Information Technologies Inc.